

PG-Seq™ A novel complete NGS solution

Customising the limit of detection for PGT-A and PGT-SR using NGS

Unlike traditional DNA analysis platforms used for Preimplantation Genetic Testing (PGT), Next Generation Sequencing (NGS) technologies provide a highly flexible and customisable workflow for the detection of segmental copy number variations (CNVs) offering variable resolution. However, greater resolution is achieved at the expense of sample throughput.

Aim – We aim to define the level of resolution achievable using the PG-Seq™ kit with different NGS sample throughput.

Materials – PG-Seq™ kit NGS files reflecting the expected total mapped read count per sample for different throughput scenarios from 1-96 samples in a single MiSeq® instrument (Illumina®) run were created *in silico* from PG-Seq™ kit NGS files using random sampling. Files with total mapped read counts of 250,000 (96 samples), 500,000 (achieved using the standard PG-Seq™ kit workflow for 48 samples), 1 million (24 samples), 3 million (8 samples), 12 million (2 samples) and 24 million (1 sample) were created. Artificial CNVs, including duplications and deletions ranging from 100 Kbp to 10 Mbp were added *in silico* across the q arm of chromosomes 1, 5, 11, 15, 17, 19, and 22 (eg Figure 1). PG-Seq™ software was used to analyse CNVs with uniform settings across all samples.

Results – As expected, the CNV resolution increased with increasing read count per sample (Figure 2).

High resolution CNV calling, for the detection of germline or *de novo* 100 Kbp duplications or deletions, required > 12 million reads, or no more than 2 samples per MiSeq® instrument run.

A standard 48 sample PG-Seq™ kit run detected 5 Mbp segmental CNVs at a frequency of > 99.5% and should theoretically be suitable for detecting *de novo* CNVs of this size.

A higher throughput of 96 samples per run detected 10 Mbp CNVs with a frequency of 100% (Figure 3). While germline segmental CNVs could be called at this low resolution, there is a risk of false positive detection of *de novo* segmental CNVs.

The location of CNVs and their proximity to other CNVs, centromeres, and telomeres affects software CNV detection. CNVs distal to the centromere were detected at a higher frequency than those located proximal to the centromere. Additionally, at lower read counts per sample, two or more nearby duplications or deletions were sometimes detected as a single CNV.

K. Warren, S. Myers, M. Jasper.
RHS Ltd, Thebarton, Australia

RHS Ltd, 40-46 West Thebarton Road, Thebarton
South Australia 5031, Australia | rhsc.com.au




Figure 1. – Chromosome 15 CNV overview showing five duplications and five deletions ranging from 100 Kbp, 200 Kbp, 500 Kbp, 1Mbp to 2 Mbp replicated in duplicate for a sample with 24,000,000 reads, equivalent to a single sample throughput NGS run.




Figure 2. – Resolution attained using the PG-Seq™ kit for different NGS sample throughput scenarios (n ≥180 per CNV size) A) 96 down to 1 samples (0.25-24 million reads) and B) highlighted section of A replicated, 96 down to 8 samples (0.25-3 million reads) per NGS run.




Figure 3. – Detection of a single 10 Mbp duplication and 10 Mbp deletion on Chromosome 15 for a sample with 250,000 reads equivalent to a 96 sample throughput NGS run.

Conclusions

In silico analysis of Next Generation Sequencing data highlights:

- technical limitations and financial implications for the detection of small segmental CNVs,
- the PG-Seq™ kit offers scalable and economical sequencing for PGT-A and PGT-SR, and
- a standard 48 sample PG-Seq™ kit NGS run can detect *de novo* CNVs as small as 5 Mbp.

These *in silico* observations should be validated using clinical samples with known CNVs, particularly if considering the detection of small *de novo* aberrations.